GUIDELINES FOR USE RIGID INTRACARDIAC SUCTION WAND S-099, S-033, S-033-90, S-022, S-077

DESCRIPTION

The RMI Rigid Intracardiac Suction Wand consists of a plastic tip attached to an angled stainless steel tube with a rigid plastic handle located on the tube. The proximal connector is sized to fit 1/4" flexible tubing. Each suction wand is individually packaged sterile and non-pyrogenic in a sealed, peel-type pouch.

INDICATIONS FOR USE: The Rigid Intracardiac Suction Wand is intended for use to remove excess fluid from the surgical field.

CONTRAINDICATIONS FOR USE: This device is not designed, sold, or intended for use other than indicated.

WARNINGS AND PRECAUTIONS

- Suction lines should be as short as possible to avoid hemolysis caused by blood pooling or mixing of blood with air.
- Sucker blood should not be reused due to its relatively high level of hemolysis and particulate debris, unless the operating surgeon deems it necessary.
- If the decision is made to reuse it, the sucker blood should be filtered to remove debris before reaching the oxygenator.
- Blood which is to be returned to the extracorporeal circulation must be adequately heparinized to prevent clotting.
- Blood should not be allowed to pool in the pericardium or the pleura as this may cause increased hemolysis.
- This device is for one-time use only and must be discarded after use.
- Do not use this device if package is damaged or has been opened.
- Do not use suction wand if it shows signs of damage.
- Federal law (U.S.A.) restricts this device to sale on or by the order of a physician.
- Only one sucker should be used on each suction line.
- The sucker lines should not be interconnected before reaching the suction pump since this would cause reduced efficiency and possible loss of suction.
- To minimize hemolysis, suction pumps should be operated at the lowest possible speed consistent with good surgical technique.
- Proper surgical procedures and techniques are necessarily the responsibility of the medical profession. These Instructions are furnished for informational purposes only. Each surgeon must evaluate the appropriate use of this device, case by case, based upon medical training, experience, and the type of surgical procedure employed.
- Dispose of used product in accordance with established hospital protocol for hazardous waste.

Edwards Lifesciences is a trademark of Edwards Lifesciences Corporation.

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Edwards Lifesciences Deutschland GmbH